

HAIYE



HY-BYH series Electronic Nasopharyngo Laryngoscope

/HY-BYH2.0, HY-BYH2.5, HY-BYH3.9, HY-BYH4.5, HY-BYH5.3/



Feature and Configuration

Product production adheres to the spirit of craftsmanship, using sufficient materials and fine workmanship, and using a smaller outer diameter to achieve a larger instrument channel;

Product structure: 3 sizes monitors are optional, 3.5 inch, 12.1 inch and 13.3 inch are touchable screen, the tip part, the operating part and the insertion part; The products are divided into 4 specifications: HY-BYH2.0, HY-BYH2.5, HY-BYH3.9, HY-BYH4.5, HY-BYH5.3



3.5 inch
touch screen
image processor



13.3 inch
touch screen
image processor



12.1 inch
non-touch
image processor

Parameter

| | |
|-----------------------|-------------------------------------------------|
| Field angle | 120°, tolerance ±10% |
| Field depth | 3mm~150mm |
| Rotation | up 0° to 130°, down 0° to 130°, tolerance ±10% |
| Working time | ≥240 mins |
| Instrument channel | 1.2mm~3.0mm (≥Ø1.0mm, HY-BYH2.0 not compatible) |
| Length of insert part | 360mm, tolerance ±5% |



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com