

HAIYE



HY-LGY

Medical Endoscopic Cold Light Source



Application

This product is a light source for endoscopes in medical institutions when using endoscopic camera systems.

Configuration

- > Strong magnetic anti-interference, stable light output
- > 30,000 hours of service life = 5 years
- > 120W high brightness, high power consumption lamp beads
- > 5800k-6500k color temperature=natural light
- > 80/92 color rendering index = high reproduction degree
- > 0-100 digital display = more intuitive
- > Knob adjustment=convenient
- > High temperature resistant fiber optic, long time corrosion resistance.

Parameter

Rated power	120W
Luminous flux	800 lm
Color temperature	3000-7000K
Brightness adjustment	0-100
Luminous area	7.7mm
Color rendering index	≥90
Maximum current	0.8A
Power	AC110-240V
Size	L336mm*W330mm*H116mm
Weight	5.5 kg
Temperature	0°-40° working environment
Bulb lifetime	30000h



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com