


# HAIYE



## HY-YPG Single-use Flexible Electronic Cystoscope



### Application

Used in conjunction with the electronic endoscope image processor produced by our company, it can be inserted into the urethra and renal pelvis through the urethra for observation, diagnosis, photography and treatment of the urethra and renal pelvis.



### Configuration

The product is a disposable handheld flexible electronic endoscope consisting of an electronic cystoscope, aa perfusion plug and an optional accessory waterproof cover.

image processor




12.1inch



13.3inch touch screen

video converter



### Parameter

Field angle	120°, tolerance ±15%
Feild depth	3~150mm
Camera pixels	400*400
Specification	HY-YPG3.9,HY-YPG4.5,HY-YPG5.3
ID of insert part	≥1.2-2.0mm ±0. 2mm
OD of insert part	≥3.9-5.4mm ±0. 3mm



**Shenzhen Haiye Medical Technology Co.,Ltd.**

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com