


HAIYE



HY-TC2 Multifuntional Airway Trolley



Application


Equipped with 12.1"HD1 electronic endoscope image processor produced by our company, it can be used with our video laryngoscope, disposable uroscope, bronchoscope and other products for clinical tracheal intubation or treatment.

It is suitable for ENT, operating room, ICU, pediatrics, emergency center, ambulance, respiratory medicine, infection department, thoracic surgery, urology and so on.

Parameter

Screen	12.1" color display
Resolution	1024×768
Field angle	89°
Function	HDMI, USB, signal , built-in storage card, support data storage and export, take photo, video,freeze frame
Compatible products	Video laryngoscope, video stylet, electronic bronchoscope, disposable ureterscope, etc
Battery	6600mAh
Rated voltage	7.2V
Charging voltage	8.4V
Battery life	≥4h
Charging time	≤5h
Battery lifetime	500 cycles ≥80% of rated capacity
Accessories	adapter, power cord, data cable, manual

Scenario



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com