

HAIYE



## HY-YM2.8

### Single-use Flexible Electronic Ureteroscope

HY-YM2.8A, HY-YM2.8B



### Application

Used in conjunction with the electronic endoscope image processor produced by our company, it can be inserted into the urethra and renal pelvis through the urethra for observation, diagnosis, photography and treatment of the urethra and renal pelvis.



### Configuration

The product is a disposable handheld flexible electronic endoscope consisting of an electronic ureteropyeloscope, a perfusion plug and an optional accessory waterproof cover.



### Parameter

Field angle	120°, tolerance ±15%
Field depth	3~50mm
Rotation	up/down 270°, tolerance ±15%
OD of insert part	2.8mm
ID of insert part	≥1.0mm
Length of insert part	700mm



**Shenzhen Haiye Medical Technology Co., Ltd.**

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: [info@haiyemedical.com](mailto:info@haiyemedical.com)

Website: [www.hyendoscope.com](http://www.hyendoscope.com)