

HAIYE



HY-ZQ series Electronic Bronchoscope

/HY-ZQ2.0, HY-ZQ2.8, HY-ZQ4.5, HY-ZQ5.8/

Usage

Used in conjunction with the electronic endoscope image processor produced by our company, it is used by medical units to provide the necessary images for observation, diagnosis, photography and treatment of the trachea and bronchi of adults and children over 3 months old. Product structure: The product is mainly composed of bending part, inserting part, operation part, conical sheath and video connection line.



Feature and Configuration

Product production adheres to the spirit of craftsmanship, using sufficient materials and fine workmanship, and using a smaller outer diameter to achieve a larger instrument channel;

Product structure: 3 sizes monitors are optional, 3.5 inch, 12.1 inch and 13.3 inch are touchable screen, the tip part, the operating part and the insertion part; The products are divided into 4 specifications: HY-ZQ2.0, HY-ZQ2.8, HY-ZQ4.5 and HY-ZQ5.8.



3.5 inch
touch screen
image processor



13.3 inch
touch screen
image processor



12.1 inch
non-touch
image processor

Parameter

Field angle	120°, tolerance ±10%
Field depth	3mm~150mm
Rotation	up 0° to 180°, down 0° to 130°, tolerance ±10%
Working time	≥240 mins
Instrument channel	1.2mm~3.0mm (≥Ø1.0mm, HY-ZQ2.0 not compatible)
Length of insert part	600mm, tolerance ±5%



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com