

HAIYE



HY-ZQ series Electronic Bronchoscope

/HY-ZQ2.0, HY-ZQ2.8, HY-ZQ4.5, HY-ZQ5.2, HY-ZQ5.8/



Feature and Configuration

Product production adheres to the spirit of craftsmanship, using sufficient materials and fine workmanship, and using a smaller outer diameter to achieve a larger instrument channel;

Product structure: 3 sizes monitors are optional, 3.5 inch, 12.1 inch and 13.3 inch are touchable screen, the tip part, the operating part and the insertion part; The products are divided into 4 specifications: HY-ZQ2.0, HY-ZQ2.8, HY-ZQ4.5, HY-ZQ5.2 and HY-ZQ5.8.



3.5 inch
touch screen
image processor



13.3 inch
touch screen
image processor



12.1 inch
non-touch
image processor

Parameter

Field angle	120°, tolerance ±10%
Field depth	3mm~50mm
Rotation	up 0° to 180°, down 0° to 130°, tolerance ±10%
Working time	≥240 mins
Instrument channel	1.2mm~3.0mm (≥Ø1.0mm, HY-ZQ2.0 not compatible)
Length of insert part	600mm, tolerance ±5%



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com