

HAIYE



HYHJ-1320 Disposable Video Laryngoscope



Application


It is provided for clinical use to provoke the patient's epiglottis to expose the glottis, and to guide medical staff to perform tracheal intubation and implement mechanical ventilation.

Suitable for anesthesia departments, ICU, pediatrics, emergency centers, ambulances, etc. in medical institutions.



Feature and Configuration

The product adopts the structural design combining stainless steel and polymer material, with light weight, high strength, full functions, long life; the product is equipped with camera of 2 million pixels, 3" color display screen, standard with three specifications of camera system, corresponding to three specifications of disposable laryngoscopes blades, can be applied to all people.



Parameter

Monitor	3" color display
Camera	Pixel: ≥ 2 million; illuminance: ≥ 3000 lux
Battery	Battery life: ≥ 6 h; charging time: ≤ 4 h; capacity: 3400mAh
Function	Take pictures, video, freeze, store and export data
Storage	-40°C~70°C ambient temperature
Humidity	$\leq 93\%$ rh



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com