

# HAIYE



## HYHJ-KC Reusable Video Laryngoscope



### Application

It is provided for clinical use to provoke the patient's epiglottis to expose the glottis, and to guide medical staff to perform tracheal intubation and implement mechanical ventilation.

Suitable for anesthesia departments, ICU, pediatrics, emergency centers, ambulances, etc. in medical institutions.



### Feature and Configuration

The host display and camera parts are connected by plug and pull mode, which is efficient, strong and beautiful. The product can be used with 6 different specifications of blades can be reused after disinfection, its specifications are Miller 0, Miller 1, Mac 2, Mac 3, Mac 4, D-blade, can be applied to all people.



3.5 inch  
(touch screen)



Miller 0



Miller 1



Mac 2



Mac 3



Mac 4



D-blade

### Parameter

|                    |   |
|--------------------|---|
| Monitor            | 3" LCD display / 3.5" color touch screen display                    |
| Camera             | pixel: ≥3million; light illuminance: ≥3000 lux                      |
| Battery            | battery life: ≥4h; charging time: ≤5h; capacity: 3400mAh            |
| Interface Function | HDMI, Type-C; Photo, video, freeze, playback, store and export data |



**Shenzhen Haiye Medical Technology Co.,Ltd.**

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: [info@haiyemedical.com](mailto:info@haiyemedical.com)

Website: [www.hyendoscope.com](http://www.hyendoscope.com)