

HAIYE



HYHJ-KP Disposable Video Laryngoscope /HYHJ-KP3C, HYHJ-KP5C/



Application

It is provided for clinical use to provoke the patient's epiglottis to expose the glottis, and to guide medical staff to perform tracheal intubation and implement mechanical ventilation.


Suitable for anesthesia departments, ICU, pediatrics, emergency centers, ambulances,



Feature and Configuration

The monitor and camera parts of the product adopt plug-and-play connection mode, which is highly efficient, sturdy and beautiful; The product is divided into HYHJ-KP3C and HYHJ-KP5C two models, respectively corresponding to Miller00, Miller0, Miller1 and Mac1, Mac2, Mac3, Mac4, Mac5, there are 8 kinds of disposable laryngoscope blades, can be applied to all people.

Single-use laryngoscope blade



Parameter

Monitor	3" LCD display / 3.5" color touch screen display
Camera	pixel:≥3million; light illuminance:≥3000 lux
Battery	battery life:≥4h; charging time:≤5h; capacity:3400mAh
Interface Function	HDMI, Type-C; Photo, video, freeze, playback, store and export data
Storage	-40°C~70°C ambient temperature
Humidity	≤93% rh



Shenzhen Haiye Medical Technology Co., Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com