


# HAIYE




## SH-X13K Endoscope Camera System (HD4K) SH-X13 Endoscope Camera System (1080P)



### Application

This product is a medical institution in the endoscopic diagnosis or treatment, need to be connected to the endoscope in advance, the endoscope observed by the human body cavity of the field of view area through the product for image acquisition, processing, storage and transmission to the monitor, and can provide endoscopic light source for the endoscope.



### Camera Parameter

Sensor	SONYCMOS
Resolution	4K/1080P
Freeze white balance	support
Digital amplification Brightness	0-8x, 1-30 steps
Endoscope	Otorhinolaryngoscopy/Laparoscopy/Cystoscopy/Hysteroscopy/Arthroscopy/Intervertebral Fiberoptic Surgery/Fiberscope
Language	English, Japanese, Chinese (Simplified), Chinese (Traditional), Korean
Working temperature	0°C~40°C

### Image System Parameter

Output interface	HDMI/DVIVGA/SDI/CVBS
Input signal	4K, 1080P
Power	AC110V~240V
USB flash disk	USB2.0/USB3.0 (Up to 32G), Video and photo
Shortcut key	Freeze/White balance/Digital zoom
Control board	Control mother board menu and shortcut keys via RX/TX

### Configuration

Standard imaging system host and camera, as well as power adapter, etc.

**Optional accessories:** trolley, medical monitor, endoscope cold light source, printer, computer and related software, etc.



**Shenzhen Haiye Medical Technology Co.,Ltd.**

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com