

HAIYE



SH-XG15K

Endoscope Camera System (HD4K)


SH-XG15

Endoscope Camera System (1080P)



Application

It is used for medical institutions to connect the endoscope to the product during endoscopic diagnosis or treatment, and the area observed by the endoscope can be captured, processed, stored and transmitted to the monitor through the product, and at the same time provide lighting for the endoscope.



Standard imaging system host and camera, as well as power adapters.

Optional accessories: trolley, medical monitor, printers, computers and related software, etc

Camera Parameter

Sensor	SONYCMOS
Resolution	4K,1080P
Freeze white balance	support
Digital amplification brightness	0-8x,1-30 steps
Endoscope mode	Otorhinoscopy/laparoscopy/cystoscopy/hysteroscopy/arthroscopy/foraminoscopy/fiberscopy
Language	English, Japanese, Chinese (Simplified), Chinese (Traditional), Korean

Imaging System Parameter

Output interface	HDMI/DVIVGA/SDI/CVBS
Input signal	4K,1080P
Power	AC110V~240V
USB flash disk	USB2.0/USB3.0 (Up to 32GB) Video and photo
Shortcut key	Freeze/White Balance/Digital Zoom
Control board	Control motherboard menu and shortcut keys via RX/TX

Light Source Parameter

Light source type	LED
Rated power	60W
Color temperature	4500-6500K
Brightness adjustment	Support
Light emitting area	3.8mm
Lifetime	3000h



Shenzhen Haiye Medical Technology Co.,Ltd.

Add: Grand Equipment Industrial Park, Shenzhen, Guangdong, China

Tel: +86-29-87872952

Email: info@haiyemedical.com

Website: www.hyendoscope.com